



Save-A-Switch

Hold-A-Callbell

Customers

Order Form

Contact Us

1.800.733.1444

SAVE-A-SWITCH T



FEATURES and BENEFITS

Patients have all chains at an arm's length.

Installation ready.

No tools required. Installs in seconds.

Low initial cost (less than repair cost)

Easy to see.

Fits all standard bed rails (round to square)

Can be used over and over (nothing to replace)

Durable 100% nylon hook and loop (can be disinfected)

Helps prevent against shock hazard.

Attaches to any chain. Comes with a 48 in. white nylon string with #10 ABS balls and a ABS connector for #6,#8 or #10 ball chain.

To prevent damage to switches, a release of 5lbs. (+/-2lbs) is part of the design.



How it works:

When the rail is lowered or the bed is moved. The two straps are pulled apart by the rail and are no longer secured together.



Home | Save-A-Switch | Hold-A-Callbell | Customers | Order Form | Contact Us | Blog
© COSCH CO 2010 | Site Design By: Business Network Solutions, LLC