

# *Bite-or-Puff*™

## Air-Operated Switch

Operate: Nurse Call, Computers, Communicators, etc.  
With Bite, Puff, or Squeeze

*For people with high-level paralysis*



*A gentle bite on the bulb  
will trigger the switch.*

*Bulb may comfortably be left in the  
mouth for an extended period of time.*

*Free 30-Day Evaluation  
Contact us for information*



*Made in U.S.A.*

### DESCRIPTION

Easily operate Nurse Call Systems, Computers & Communicators (with single-switch capability), and other low power devices with ***Bite-or-Puff***. The pressure required to operate it is very small and is provided by a gentle bite, puff, or squeeze on a bulb connected to the input tube.

The bulb can be a sterile *Bite-or-Puff* bulb (20 included), an infant ear syringe, or even a child's squeeze toy. This versatile general-purpose switch can also be operated as a puff switch by cutting off the end of the sterile bite bulb with sharp scissors. Bite frequently has a significant advantage over puff, since (1) Jaw muscles are powerful, (2) Bite does not need lip seal around the bulb, and (3) It is easier to keep clean. Often, people who can't use any other switch can use ***Bite-or-Puff***.

### SPECIFICATIONS

**Maximum Pressure:** 25 psig

**Trip Pressure:** Preset to 1 inch of water (.04psi)

**Switch:** Momentary Contact, gold plated contacts

**Max. Switch Rating:** 25 Volts, 25 mA AC/DC

**Connector & Cable, (Please Specify):**

- **3.5 mm "Mini Phone Plug"** with 3 foot cord for computer or communicator use- **or**
- **1/4 inch "Standard Phone Plug"** with 10 foot cord for Nurse Call use.

Other connectors & cords available on special order.

**Air Tube:** Silicone Rubber tubing, 0.1" ID by 46".

**Bite Bulbs:** Sterile Pack, 20 each, 0.4" dia. x 3.5"

**Warranty:** 1 year, Parts and Workmanship.

**MED LABS, INC.** 7334 G Hollister Ave., Goleta, California, 93117, USA  
Voice or FAX: 805-968-2486; Toll Free: 800-968-2486 [www.medlabsinc.com](http://www.medlabsinc.com)

© 2015, Med Labs, Inc.