



The **Presscall®** is used by patients with sufficient motor skills to either press or lean on the bulb to activate a call. It has an oval shaped, soft vinyl bulb design that prevents it from rolling off the patients bed. The product is **oxygen safe** and available in various colors and lengths.

Technical Specifications:

- "**Non Powered**" Pneumatically Operated Switch
- Operating voltage of up to 24 volts AC or DC at .33 amps.
- Standard contact type: Normally open.
- Standard size: Six foot tube 5/16" X 3/16" diameter.
- Standard color: Grey.