



# Health Care Call System

## Callcords and Dummy Plugs Cat. 8833/8835/8836



### Technical Data

#### Pendant and Cord:

Grey molded thermoplastic pendant with grey vinyl cord

#### Cord Length:

8833: 6 ft. (1.8m)

8835: 10 ft. (3m)

#### Controls:

Call Button: Momentary contact

#### Connector:

1/4 in. (0.64cm) plug

#### Bed Clip:

Removable white nylon strap with spring-type, plated, metal clip.

#### Certification:

UL 1069 Listed

C-UL – Canadian Standard

C22.2 No. 125

### Description

Jeron's Cat. 8833/8835 Callcords and Cat. 8836 Dummy Plugs are used with patient stations in EC-300 and Pro-Alert™ tone-visual nurse call systems.

### Features

- Heavy-duty construction
- Strain relief at both ends
- Spring clip for attachment to bedding

### Operation

Pressing the button places a call on the system. A cordout call will be placed on the system if the callcord is removed from the patient station jack. The cordout call is cancelled by inserting the callcord plug back into the patient station jack. A Cat. 8836 Dummy Plug may be inserted into the patient station jack to prevent cordout calls from being placed when a callcord is not connected at the patient station.

### Architects/ Engineers Specifications

Callcords shall be Jeron Cat. 8833 (6 ft.) or 8835 (10 ft.), made of grey thermoplastic with momentary contact Call pushbutton, vinyl insulated cord, metal bed clip, strain relief, and 1/4" plug.

### Associated Jeron Equipment

- Cat. 6121/23 Patient Station

PRODUCT SUBJECT TO CHANGE WITHOUT NOTIFICATION