



# Model No. 4T-1B

## Single Patient Station

### DESCRIPTION:

The 4T-1B is a patient bedside station that is intended to be wall mounted at the headwall above a patient bed. The station includes one call cord receptacle which will accept a variety of 1/4" jack momentary call cord sets. The station also includes an individual call indicator LED to show the presence of an active call.

When the call cord is pressed the call will illuminate in three locations simultaneously. The stations' call assurance indicator, the rooms' corridor dome light, and the enunciator panel at the nurses' station. In addition, a slow repeating call tone will sound at the nurses' station and any associated duty stations.

Removal of a call cord set will automatically initiate a call. The call cord or a dummy plug must be replaced before the call can be reset.

The station can only be reset at the point of origin unless an in-room remote reset station has been added for your convenience.

### OPTIONAL ACCESSORIES:

CS-1 Momentary Cord set

DP-1 Dummy Plug

4T-RR Remote Reset station

B09050-04 Connector pigtail



### FEATURES:

- Durable faceplate
- Individual bed indicator
- Plug-in wiring header for easy removal
- Scratch resistant Lexan® overlay
- Cord disengaged safety feature

### SPECIFICATIONS:

**Circuit Board:** PCB mounted header for use with insulation displacement connectors. Note: Connectors sold separately.

**Back Box:** Single gang 2-1/2" deep. Steel City #CD or equal.

**Dimensions:** 4.50"H x 2.75"W x 2.00"D

**Faceplate:** Constructed of UL94V-0 ABS flame resistant plastic.

**Switch:** Momentary activation, mechanical life >1 million operations.

**Hardware:** Includes stainless steel mounting screws.