



Call Crest Technical Support for assistance at: **1-800-328-8908**

www.cresthealthcare.com

©2010 Crest Healthcare Supply



Patient Station, Crest Replacement for Dukane

4A1220C

4A1225C (With Speaker)

Function Table	
COLOR VARIATIONS	
DUKANE	CREST
Purple with Yellow stripe	Purple
Blue with Yellow stripe	Blue with White stripe
White with Red stripe	Red with White stripe
Black	Black
White	White
Yellow	Yellow
Red	Red

Typical Installation

Please refer to your original Dukane installation manual 403-56 or 403-57 for external connections.

Specifications

ELECTRICAL

Current and Voltage Rating: 24VDC -80MA

MECHANICAL

Mounting: Fits standard single gang electrical back box

SIZE MEASUREMENTS

Plate: 4.5" x 2.75" (1 Gang)